

C-Therm TCi™ Thermal Conductivity Analyzer

*Introduction to the
Modified Transient Plane Source
(MTPS) Technique*



Agenda

- Who we are?
- What we do?
 - Thermal Conductivity (k)***
 - ✓ Introducing to the TCi™
 - ✓ Principles of Operation
 - ✓ Capabilities
- Applications
 - ✓ Metal Hydrides
 - ✓ Explosives
 - ✓ Polymers & Nanomaterials
 - ✓ Electronic Adhesives
 - ✓ Carbon Nanotubes (CNTs)
 - ✓ Asphalt / Concrete
- Comparison of Methods
- Q&A



C-THERM
TECHNOLOGIES ^{Ltd.}

WHO WE ARE?

Who is C-Therm?

C-THERM
TECHNOLOGIES ^{Ltd.}

(Formerly Mathis Instruments Ltd.)

Non-destructive thermal sensor technology solutions for R&D, production, and QC applications, delivering fast, accurate measurement of **thermal conductivity** and **effusivity** in seconds with virtually unlimited sample size.





C-THERM
TECHNOLOGIES^{Ltd.}

WHAT WE DO


C-Therm Product Lines

THERMAL CHARACTERIZATION

C-Therm TCi™ Thermal Conductivity Analyzer

Clients include:

- NRC
- Philip Morris
- US Army
- US Navy
- Kodak
- Texas Tech University
- P&G
- Henkel Technologies
- Dupont
- Atomic Weapons Est. (AWE)



PHARMACEUTICAL APPLICATIONS

C-Therm ESP™ Effusivity Sensor Package Process Analytical Technology (PAT)

Clients include:

- Patheon
- Astra Zeneca
- Wyeth
- Biovail
- BMS
- USP

Overview of TCI™

Fast

0 to 100 W/mK in 5 seconds

Wide Temperature Range

-50° to 200°C

Easy-to-Use

No calibration required

No Sample Preparation

Unlimited sample sizes

Non-Destructive

Leaves sample intact

Versatile

Tests solids, liquids, powders & pastes

Highly Flexible

Designed for lab, QC & at-line testing

Accurate

Better than 5%

Precision

Better than 1%



What does it Measure?

The C-Therm TCi measures two thermal properties primarily:

$$\text{Thermal Conductivity} = (W / m \cdot K)$$

and

$$\text{Effusivity} = \sqrt{k\rho c_p}$$

Where : k = Thermal Conductivity (W/m • K)

ρ = Density (kg/m³)

c_p = Heat Capacity (J/kg • K)

It also indirectly measures (calculated) Thermal Diffusivity and Heat Capacity and has user input capabilities to determine Density

How the sensors work?

Wood feels warm



Metal feels cold



Heat always flows from a hot object to a cold object.

Wood is not a good conductor of heat, so it is **slow** to absorb the heat.

Metal has higher “**thermal conductivity**” so the heat from your hand flows into the metal **quickly** - creating the sensation of it being cold.

C-Therm sensors work like your hand, by **rapidly** determining the **rate** of heat flow from one material to another. Like your hands, our sensors **supply** the heat source *and* **detect** the heat flow. They also have no **sample size** issues, and do not destroy the sample being tested.



How the Technology Works


Sample material can be solid, liquid, powder or paste.

A known current is applied to the sensor's heating element, providing a small amount of heat.

The heat provided results in a rise in temperature at the interface between the sensor and the sample – typically less than 2°C.

This temperature rise at the interface induces a change in the voltage drop of the sensor element.

The rate of increase in the sensor voltage is used to determine the thermo-physical properties of the sample material.




How the sensors work?

The thermal conductivity of the sample material is inversely proportional to the rate of increase in sensor voltage. The change in voltage drop correlates with an increase in temperature at the sensor interface.

The more thermally insulative the material is – the steeper the voltage rise.

Results are displayed on the system's laptop computer in real time.




First 0.3 Seconds:
Addressing Contact
Resistance, Non-
Linear

0.3 – 0.8 Seconds:
Within Sample, Linear

No Interpretation Required

Results are presented in tabular format displaying:



TCI - [Test Results: 146 C-THERM ACH]

TCI File Manage View Reports Security Tools Windows Help

New Test Open Test Manage Diagnostics Reference Test Log Off

Test Results Filter

Test Details Test Notes Show Samples Show Settings Show Derived Values

All Tests

- Show Tests
- by Project (3 found)
- Default
- Trelleborg
- UNB
- Show Tests
- by Material (17 found)
- by Instrument (3 found)
- by TestMethod (3 found)

#	Sensor	Valid	Effusivity (W/s/m ² K)	k (W/mK)	1/m	R2	Stat	VO (mV)	VMax (mV)	DeltaV (mV)	Ambient (°C)	T0 (°C)	DeltaT (°C)
1	T067	<input checked="" type="checkbox"/>	203.0	0.088	61.64	0.9978	17:44:09	2419.3	2428.5	9.2	29.5	28.8	1.2
2		<input checked="" type="checkbox"/>	194.9	0.086	61.16	0.9979	17:45:15	2417.6	2426.9	9.3	29.7	28.7	1.2
3		<input checked="" type="checkbox"/>	193.4	0.085	61.09	0.9977	17:46:21	2417.6	2426.9	9.3	29.8	28.8	1.2
4		<input checked="" type="checkbox"/>	194.6	0.086	61.09	0.9977	17:47:26	2418.5	2427.8	9.3	29.3	28.7	1.2
5		<input checked="" type="checkbox"/>	194.7	0.086	61.15	0.9977	17:48:32	2417.4	2426.7	9.3	29.8	28.7	1.2
6		<input checked="" type="checkbox"/>	194.9	0.086	61.17	0.9979	17:49:37	2417.5	2426.8	9.3	29.8	28.7	1.2

Thermal Effusivity

Thermal Conductivity

Temp.

Easily exportable to EXCEL™

Fast


- No Calibration or Sample Preparation
- Instrument comes fully calibrated from the factory
- No sample preparation required for testing
- Rapid Testing
- 0 - 100 W/mK in only 0.8 to 5 seconds
- Accelerates your research or QC processes dramatically
- Highly accurate, repeatable testing
- Add a Second Sensor, theTCi can be equipped with two sensors to double your throughput



Flexible

Testing Solids


- Breakthrough in the characterization of critical performance attributes of ceramics used for aerospace applications
- Simplicity of sample format is key advantage
- Eliminates technician time required for sample preparation
- Sample size flexibility evaluation of actual product formats
- Eliminates the need to mock-up samples



Flexible

Testing Liquids


- **Helping manufacturers improve heat transfer properties of advanced nano-filled liquids**
- **Allows wide range of operating temperatures for the testing of engineered liquids**
- **Low amount of heat introduced during testing minimizes convective errors typical of liquid testing**
- **Can measure liquids through a bag – or directly**



Flexible

Testing Powders


- Only instrument engineered for evaluating the thermal conductivity of powders
- Used to safely test the stability, degradation, and shelf life of explosives
- Sample volumes can be as small as 3/8 tsp
- Minimal samples are critical to pharma customers, where active ingredients are expensive and available in low quantities
- Migratable to manufacturing as a cost-effective way to monitor powder processes



Flexible

Testing Pastes

- Providing insights into all materials that contribute to the overall thermal budget for microprocessors
- Thermal interface pastes are a prime application for the TCI
- Testing allows different amount of heat penetration – resulting in a variable scale of scrutiny to ensure homogenous distribution of vital filler components



Flexible

Operate the TCi sensor in various environmental enclosures:


Thermal Chambers



Glove Boxes



Easy




There isn't a faster, easier way to generate thermal conductivity results on your sample materials!!!!

Step 1. Place sample on sensor.

Step 2. Press Start.

Step 3. Wait less than 5 seconds for your first result.



C-THERM TciTM
Thermal Conductivity Analyzer

APPLICATIONS

Applications

IMPORTANT NOTE: with applicability in a wide and diverse range of industries, the following section touches only on a select few highlighted applications.

For more information on the comprehensive range of applications for the TCI please visit www.ctherm.com and refer to our vast technical library.



Application: Fuel Cell (Hydrogen Storage)

Background



Pic: Dr. Shane Beattie of the University of New Brunswick working with the TCi

In the development and characterization of metal hydrides for hydrogen absorption, the **thermal conductivity** of the material is an important material performance property. The thermal conductivity of the solid material is **important to know how long it will take to get the pellets to an equilibrium temperature to understand when they will begin releasing hydrogen.** Furthermore, it is important to know the thermal conductivity of these pressed powders while **engineering a storage device to use the powders commercially.**

The thermal conductivity of the material changes substantially from a powder to a compressed solid necessitating the direct measurement of both material formats.

Application: Hydrogen Storage

Special Setup

Due to the highly reactive nature of the materials the sample are tested in an argon or nitrogen purged glove box.


The TCi sensor cable is fed through one of the sealed ports – enabling the sensor to test both the powder and solid pressed pellets inside of the glove box.



Application: Hydrogen Storage

Results

Thermal Conductivity of Metal Hydrides



Application: Explosives & Energetics

■ Application


Safety constraints have traditionally presented researchers with challenges in testing the thermal conductivity of energetic powdered materials. Minimal volumes of energetic powders can now be tested with the TCI in employing the modified transient plane source technique.

■ Profile Clients

- ✓ US Army & Navy (2X Each)
- ✓ NRC (Natural Resources Canada)
- ✓ Australian Defense Department

■ Tools & References

- ✓ Harris & Sorensen. *Testing Minimal Volumes of Energetic Powders*
- ✓ Singh (et al). *Physical Characterization, Accelerating Rate Calorimetry and Thermal Conductivity of Ammonium Nitrate Emulsion*
- ✓ **NEW:** *Archived Webinar Recorded by NRC*



Frank Kamenetskii Equation

$$T_c = \left[\frac{R}{E} \ln \left(\frac{d^2 \rho H Z E}{T_c^2 \lambda \delta R} \right) \right]^{-1}$$

Where,

- d - diameter
- H,Z,E - kinetic parameters
- δ - shape factor
- λ - thermal conductivity
- R - gas constant
- ρ - density

Background

- **Critical Temperature using the Frank-Kamenetskii equation**
 - ✓ Kinetics of decomposition, heat capacity and thermal conductivity
- **Thermal conductivity via MDSC is not possible for powders/liquids**
- **Need for rapid, safe measurement**
 - ✓ Minimize volume requirement
 - ✓ Eliminate ESD potential
 - ✓ Limit opportunities for powder ingress

Accuracy

Table A1: Accuracy Test Results vs. Certified Values

Sample	Measured k (W/mK)	Avg. of the Measured k (W/mK)	RSD (%)	Real k (W/mK)	Accuracy (%)
PDMS	0.162194	0.162865	0.24%	0.159	2.43%
	0.162951				
	0.162945				
	0.163008				
	0.163226				

Reference Material: PDMS (DiMethyl PolySiloxane Silicone Fluid)

Sample Material

Ammonium perchlorate (NH_4ClO_4) was measured with a bulk particle size of 200μ



Ammonium perchlorate: used commonly as an oxidizer in solid rockets.

Testing Method




Pictures are example only and do not show ground wires.

Results

Test	Thermal Conductivity (W/m·k)	Average (W/m·k), N=4	RSD ¹ (%), N=4	Average (W/m·k), N=3	RSD (%), N=3
1	0.127086	0.1276	0.31%	0.1288	0.75%
	0.127641				
	0.127866				
	0.127969				
2	0.130499	0.1293	0.64%		
	0.128595				
	0.128970				
	0.129272				
3	0.129295	0.1293	0.38%		
	0.129967				
	0.128839				
	0.129063				

Results



Conclusion


- The thermal conductivity of a volume of approximately $\frac{3}{8}$ teaspoon of ammonium perchlorate was measured with an accuracy of better than 2.5%.
- Further studies are recommended to investigate the relationship between particle size, moisture content, and packing density on the thermal conductivity.

Application: Carbon NanoTubes (CNTs)

Background

Carbon NanoTubes (CNTs) have novel properties that make them potentially useful in many applications in nanotechnology, electronics, optics and other fields of materials science, as well as potential uses in architectural fields. **They exhibit extraordinary strength and unique electrical properties, and are extremely efficient conductors of heat.** CNTs tend to aggregate and cluster

however,



SEM images of pristine SWCNT/epoxy and shortened SWCNT/epoxy composites.
(a) Dispersion of pristine SWCNTs in epoxy resin. (b) Dispersion of shortened SWCNTs in epoxy resin.


resulting in poor dispersion and inferior properties.

Thermal conductivity is an important performance attribute as researchers

assess the impact of mechanically altering the single walled carbon nanotubes (SWCNTs) to produce short and open-ended fullerene pipes for better dispersion used in polymer composites.

Application: Carbon NanoTubes (CNTs)

Results



The University of Texas Tech found with their TCI that shortening the CNTs did not adversely impact the thermal conductivity. This was important as shortening the CNTs considerably improved the dispersion in the composites.

Application: Carbon NanoTubes (CNTs)

■ Profile Clients

- ✓ Texas Tech University



■ Tools & References


- ✓ Wang, Liang, Wang & Zhang. *Dispersion and thermal conductivity of carbon nanotube composites*
 - LINK: <http://dx.doi.org/10.1016/j.carbon.2008.08.024>

Application: Asphalt and Concrete



Application: Asphalt and Concrete

Sample



Application: Asphalt and Concrete

Specimen	Bulk density, kg/m^3	C-Therm TCI Modified Transient Plane Source (MTPS)		Guarded Hot Plate (GHP)		Traditional Transient Plane Source (TPS)	
		k, $\text{W/m}\cdot\text{K}$	RSD*, %	k, $\text{W/m}\cdot\text{K}$	RSD, %	k, $\text{W/m}\cdot\text{K}$	RSD**, %
LW1	1686	0.841	2.62	0.845	n/a	-	-
LW2	1587	0.786	3.24	0.785	n/a	-	-
NW1	2585	2.124	3.12	-	-	2.008	1.12
NW2	2459	2.206	2.07	-	-	2.187	1.16



*RSD of five measurements across six different locations.

**RSD of five measurements at the same location.

Application: Food Products

Conclusion:

- Results show a good correlation
- Sample results for the light weight asphalt (LW1 and LW2) and the normal weight (NW2) provide a better than 1% agreement.
- Sample results for NW1 differed by approximately 6% which is likely due to positional bias of where the TPS sensor was located in relation to the multiple sample locations tested with the MTPS sensor.
- Advantage to the TCi is the ability to look at the relative homogeneity of the sample in understanding how this impacts the localized thermal conductivity performance.



C-THERM TCI™
Thermal Conductivity Analyzer

HOW DO WE COMPARE

Competitive Positioning


Superior to Other Test Methods:



Mathis TCI
(Modified Transient
Plane Source)



**Traditional Guarded
Hot Plate**



**Transient Plane
Source**



**LaserFlash
Diffusivity**

Sample Configuration

Minimum	0.67" diameter (17mm)	6" x 6" (150 x 150mm)	Two Identical Samples 1" x 1" (25 x 25mm)	0.5" diameter (12.4mm) 0.004" thick (1mm)
Maximum	Unlimited	24" x 24" (600 x 600mm)	Two Identical Samples Unlimited	0.5" diameter (12.4mm) 0.004" thick (1mm)
Material Testing Capabilities	Solids, Liquids, Powders, Pastes	Solids	Solids, Liquids	Solids

Competitive Positioning


Superior to Other Test Methods:



Mathis TCI
(Modified Transient
Plane Source)



**Traditional Guarded
Hot Plate**



**Transient Plane
Source**



**LaserFlash
Diffusivity**

Speed & Flexibility

Sample Preparation	None Required	Extensive	Some	Extensive
Testing Time	Seconds	Hours	Minutes	Seconds*
Training Time	Minimal	Moderate	Significant	Extensive
Non-Destructive	Yes	No	No	No
Integrated, Downloadable Test Results Database	Yes	No	No	No

Competitive Positioning


Superior to Other Test Methods:



Mathis TCI
(Modified Transient
Plane Source)



**Traditional Guarded
Hot Plate**



**Transient Plane
Source**



**LaserFlash
Diffusivity**

Range

k-Range (W/mK)	0 – 100	0 – 2	0 – 100 (100 – 500 requires C_p)	0 – 500
Temperature Range (°F) (°C)	-58° to 392°F -50° to 200°C	-4° to 392°F -20° to 200°C	-148° to 2552°F -100° to 1400°C	-148° to 3627°F -100° to 2000°C
Pricing	\$	\$ \$	\$ \$	\$ \$ \$

Comparison to Traditional Techniques

Material	TCi Result (W/mK)	Reference Value (W/mK)	Comparative Technique	Source	% Difference
PYREX	1.17	1.135	Guarded Hot Plate (GHP)	BCR Certification Report for a Pyrex Reference Sample	3.04%
PDMS	0.162	0.159	Hot Wire (HW)	Bashar, University of New Brunswick	1.87%
LAF	0.0880	0.086	Transient Plane Source (TPS)	Hume, 3 rd Party Contract Lab	2.30%
Water -gel	0.626	0.627	Hot Wire (HW) & Modified Hot Wire (MHW)	Picot & Bashar, University of New Brunswick	0.16%


Presented by Dr. Robert Bateman of the University of New Brunswick at ITCC (International Thermal Conductivity Conference) 2009, Pittsburgh PA

Client Testimonial



Website: www.ctherm.com

- Resources
- Technology
- Video Demonstration
- Client Testimonials
- News



Upcoming Webinars

Date

Topic

Thur., April 7th

TCi Application Highlight: Rock On! (Geology)

Tue., April 12th

TCi Application Highlight: Hit the Road (Asphalt & Concrete)



Questions for C-Therm?



Adam Harris
Managing Director
C-Therm Technologies Ltd.
(Formerly Mathis Instruments Ltd.)

Email: aharris@ctherm.com

Toll-Free: 1-877-827-7623
 (North America)

Direct: (506) 461-7203

THANK YOU!