

A [printed circuit board](#) (PCB) is defined as an item which, “mechanically supports and electrically connects [electronic components](#) using [conductive](#) tracks, pads and other features [etched](#) from copper sheets [laminated](#) onto a non-conductive [substrate](#).” (definition courtesy of [Wikipedia](#)). PCBs may be single- or double-sided, and may have multiple layers. Some have components embedded to allow for complex and advanced circuitry (Figure 1).



*Figure 1. An advanced PCB with embedded components visible ([Source](#)).*



Figure 2. A PCB with a visible failure, caused by overheating. Such a failure can be due to poor thermal management ([Source](#)).

## **Importance of Thermal Conductivity with Printed**

## Circuit Boards

[Heat management](#) is crucial for PCB performance, reliability, and longevity. Inadequate heat management may lead to [delamination](#), damage, or [device failure](#) (Figure 2). [Thermal conductivity](#) plays a vital role in heat management, and thus it is a key parameter for PCB design. The [C-Therm Trident thermal conductivity platform](#) is a useful tool in obtaining [rapid, precise, and accurate](#) measurements of PCB component thermal conductivity.

Here, we detail how the [C-Therm Technologies TCi sensor](#) may be used to aid in PCB heat management design. This application note comes with a special foreword by [Doug Brooks](#), a PCB design expert, regular contributor to industry publications such as [PCB Design 007](#), and owner of [UltraCAD Design, Inc](#):


For the last 25 years I have owned a printed circuit board (PCB) design service bureau. The “hot” topic in PCB design is high-speed [signal integrity](#). But on another front, PCB designers may be interested in how hot (literally) an individual [PCB trace](#) becomes. Trace temperature is directly related to [reliability](#). In the extreme, a trace that is too hot can melt the solder or cause a board to [delaminate](#). But generally we want trace temperatures to be a lot lower than that. For very high reliability applications (e.g. manned space, medical, etc.) we may want to design very conservatively. For consumer products we can be a little more aggressive. For applications in a hot desert (think war time) we may want to know how much heat we have to dissipate through some external means. The trace reaches a stable temperature when the [heating of the trace equals the cooling of the trace](#). The heating of the trace is [caused by the  \$I^2R\$  \(power\) drop](#) across the trace. The cooling of the trace is primarily the result of [conduction](#) through the [dielectric](#) (board material), and secondarily through [convection](#) and [radiation](#). It is only in the last 10 years that the industry has recognized the importance of the dielectric in the trace cooling process. The important material property in trace heating is the [resistivity](#) of the trace material (typically copper foil or plating.) Although the *actual* resistivity of a trace is subject to some discussion in the industry, most estimates are that it is between that of pure copper (1.7  $\mu\text{Ohm-cm}$ ), and about 2.1  $\mu\text{Ohm-cm}$ . The important material property in trace cooling is the [thermal conductivity](#) (W/mK) in the x,y plane and in the z axis. This can vary significantly between material offerings and even between manufacturers for the same material specification. Moreover, not all manufacturers provide a thermal conductivity specification, particularly in the z-axis. The C-Therm Thermal Conductivity Analyzer can be an important tool in measuring the thermal conductivity of a board material or in verifying that a manufacturer’s product meets specification. For more information about the relationships between PCB

trace current/temperature see the series of five papers posted at [www.ultracad.com](http://www.ultracad.com) (with at least two more coming in future months.)

Douglas Brooks, President UltraCAD Design, Inc. [www.ultracad.com](http://www.ultracad.com)

## Case Highlight: Testing the Thermal Conductivity of an FR4 PCB Dielectric Material

Samples of an FR4-type PCB dielectric material were obtained as pale yellow rectangular prisms and used as provided. Thermal conductivity was measured using a C-Therm Modified Transient Plane Source (MTPS) technique.



In-plane thermal conductivity of an FR4-type PCB dielectric material

The MTPS technique employs a patented Guard Ring technology which results in a one-directional heat flow and, thus, thermal conductivity measurement. This enables easy measurement of anisotropic materials by positioning the sensor along different surfaces of a material to obtain direct measurement of the thermal conductivity anisotropy.

Through-plane thermal conductivity data was obtained for three samples with three tests of fifteen measurements each. The mean of the three samples was taken and reported as the through-plane thermal conductivity. In all tests, the relative standard deviation was <1%, and across the three samples, the relative standard error was 0.7%. The in-plane thermal conductivity of the bulk material was measured with three tests of five measurements each. Relative standard deviation was 0.6% across three tests for in-plane thermal conductivity.

This work illustrates the importance of obtaining thermal conductivity data along all relevant material axes in the case of materials with anisotropic thermal conductivity - in this case, measuring only the through-plane value would lead to a significant underestimation of the effective thermal conductivity of the PCB dielectric material.